

MANAGING CUSTOMIZED SYSTEMS IN A HIGHLY STANDARDIZED BIOPHARMA INDUSTRY

SUS BINGO



FAT	Manifold	Lastenheft	FMEA	ETD
Bag agnostisch	FS	Plate Freeze and Thaw units	RoSS.Fill Bag Unit	QM Documents
CAD-Zeichnung	Steril	SAT	Kick-off	URS
Kostenanalyse Entwicklungs-kosten	Product Life Cycle	ISO 9001	Single Use Support	Biopharma
Pflichtenheft	RoSS	Project Management	Design Dokumente	Project Plan

COMPANY HIGHLIGHTS





SITES IN AT & USA
DISTRIBUTORS IN AMS & AS

LOCATIONS




300+

PLATFORMS SHIPPED



100 Mio.

EUR TURNOVER 2022



200.000+

Ross® Shells Shipped



130+
EMPLOYEES




25+

COUNTRIES




WE CLOSE THE GAP



ADVANCED FLUID MANAGEMENT








CONSUMABLES

Discover our single-use consumables for optimal fluid management and product protection during filling, freezing and transportation.

PROTECTION



RoSS® Shell

Protection for any singleuse bag to ensure 0% product loss

PRIMARY PACKAGING



IRIS Bioprocess Containers

Reliable primary packaging

SINGLE-USE ASSEMBLY



IRIS Single-Use Assemblies

Fluid transfer with customized single-use assemblies

COLD CHAIN SHIPPING



RoSS.SHIP

Cold chain shipping of frozen drug substances

ASPETIC SAMPLER



IRIS.SMPL

Sampling in biopharmaceutical aseptic processes



PLATFORM SYSTEMS

Single Use Support's platform solutions for a completely secure and fast logistic process.

FILL & FILTRATION



RoSS.FILL

Modular and automated aseptic filling platform for single-use systems.

FREEZE & THAW



RoSS.pFTU

Controlled automated freezing and thawing process done safely and fast

ULTRA COLD STORAGE



RoSS.FRDG

Reliable and safe ultra-cold storage

CRYOGENIC FREEZER



RoSS.LN2F

Cryogenic controlled rate freezer for temperatures down to -180°C/-292°F

HOMOGENIZATION



RoSS.PADL

Scalable cooling and massaging platform to avoide agglomerations

HOMOGENIZATION



Ross.shak

Reliable shaking and thawing of biopharmaceuticals

COMPARISON

Static vs. Plate Freezer 10L bag

Static Freezer

"Slow freezing"

Dark blue region at the middle bottom visible



Plate Freezer

"Fast freezing"

More uniform coloring



The 10L bag in the plate freezer shows a significant homogeneous coloring image as the 10L bag which was frozen in the static freezer.

THE BIOPHARMA INDUSTRY TRUSTS IN US




















SINGLE USE SUPPORT'S END-TO-END



Development under GAMP5 regulation

Controlled & documented product life cycle for GMP customers



Produkt Management

- >Konstruktion / Technische Zeichnung
- >Tech. Abklärung mit Lieferanten
- > Erstmusterprüfungen
- >Pflegen von Stücklisten
- > Prototypenbau
- > Versuche
- >Kostenrechnungen/Schätzungen
- > Projektplan / Projekt Dokumentation

Produkt Management



